

The Complete Solution for Clean Room Aerosol-Based Disinfection

MINNCARE DRY FOG SYSTEM

Current Methods of Room Disinfection

- ⇒ Surface Wiping
- ⇒ Manual Spraying (with Spray Bottles)
- ⇒ Heating Process (Vaporization)
- ⇒ Cold Process - Wet fogging
 - **Dry fogging**

✘ Disinfection is NOT a Cleaning Process

Currently Used Chemical Products


- ⇒ Alcohol based
- ⇒ Formaldehyde based
- ⇒ Gluteraldehyde based
- ⇒ Quaternary Ammonium Compounds
- ⇒ Cocktails of the above

Current Problems

- ⇒ Efficiency (more or less broad spectrum biocide)
- ⇒ Adaptation of the Micro-Organisms
- ⇒ Toxicity of the Chemicals Used
- ⇒ Long Contact Time
- ⇒ Very Long Time for Venting
- ⇒ Neutralization Needed
- ⇒ Corrosion
- ⇒ Variability of Chemical Application
- ⇒ Residuals



Minnicare® Cold Sterilant



⇒ An alternative solution for room

disinfection in Pharmaceutical,
Biotech, Medical Facilities,...

Benefits of Minncare®

- ⇒ Superior Biocidal Activity
- ⇒ All Components are Ultra-pure, Pharmaceutical Quality
- ⇒ US EPA (Environment Protection Agency) Registered Sterilant
- ⇒ No Heavy Metal Trace Contamination / Stabilizers.
- ⇒ Biodegradable – Decomposes to Acetic Acid, Water and Oxygen
- ⇒ No Toxic Aldehyde Type Vapors, Easy to Vent
- ⇒ Validated Residual Vapor Detection System

Minncare® Applications

- ⇒ Water Systems Disinfection
Tanks, Piping, Resins, Filters, RO Membranes,...
- ⇒ Surface Disinfection (Wiping and/or Manual Spraying)
- ⇒ Fogging



Minncare® Biocidal Activity

Minncare® is a potent antimicrobial agent and is effective against a broad spectrum of microorganisms, including:

- ⇒ Bacteria
- ⇒ Yeast and Molds
- ⇒ Mycobacteria
- ⇒ Bacteria Spores
- ⇒ Viruses

Activities of the most important biocides (Guyader, 1996)							
Biocides	Bacteria		Mycobacteria	Spores	Moulds	Yeasts	Virus
	Gram -	Gram +					
Peracetic acid	+++	+++		++	++	++	++
Alcohols	++	++		0	++	++	+
Alcohol (70?)	++	++	0	+	+	++	+
Glutaraldehyde	+++	+++	++	+	+++	++	++
Quat Ammonium	+++	+*	0	0	+	+	+
Chlorine	+++	+++	++	++	++	++	++
Hydrogen Peroxyde	+++	+++		+	+	+	0
Iodine	+++	+++	++	++	++	++	++

* Not active on *Pseudomonas*

Minncare® Mechanism of Action

- ⇒ Disrupts Sulfhydryl (-SH) and Sulfur (S-S) bonds in proteins and enzymes. Important components in cells and membranes are broken by oxidative disruption.
- ⇒ Impede cellular activity by disrupting chemosmotic function of lipoprotein cytoplasmic membrane transport through rupture or dislocation of cell walls

- ⇒ Denature the properties of protein components by altering the nucleic acid structure of organisms.
- ⇒ Damage vegetative cells by oxidation with hydroxy radicals.
- ⇒ Produces organic radicals that act as reducing agents for spores.

“NO CHANCE” For Micro-Organisms to Build a Resistance

Minnicare® In Vitro Aqueous Testing

Species	Minnicare® Concentration .5% Count per ml	Time for 100% kill in minutes
Bac. subtilis	6×10^5	2.5
Bac: stearothermophilus	6×10^6	2.5
Bac. subtilis NCTC 3610	2.4×10^9	5.0
Bac. mesentericus	1.6×10^9	5.0
Clostr. perfringens	1×10^7	10.0
Clostr. tyrobutyricum	1×10^7	5.0
Sacchar. cereisiae	6×10^7	0.5
Cand. mycoderma	1.4×10^8	0.5
Hansenula anomala	6.4×10^8	0.5
Pichia membronaefaciens	4.8×10^8	0.5
Pen. cameronense	1.7×10^8	2.5
Mucor plumbeus	3×10^6	2.5
Geotrichum candidum	2×10^7	0.5
Byssoschlamys nivea	6×10^7	0.5
Staph. aureus	2.6×10^9	0.5
Strept. faecalis	4.6×10^9	0.5
Kleb. aerogenes	2.3×10^9	0.5
Ps. fluorescens	4.6×10^9	0.5
Ps. aeruginosa	2×10^9	0.5
Salm. thyphimurium	2.8×10^9	0.5
Coryneb. rubrum	1×10^7	1.0
Leuconostoc spec.	5.3×10^9	0.5
Lactob. brevis	1.8×10^9	0.5

Minnicare® AOAC Sporicidal Test

A total of 720 tests were performed on three lots of the test chemical. Results of the test are shown on the chart below.

Lot/Organism	Carrier type	Positives/Total
1 Bacillus subtilis	loop	0/60
1 Bacillus subtilis	cylinder	0/60
1 Clostridium sporogenes	loop	0/60
1 Clostridium sporogenes	cylinder	0/60
2 Bacillus subtilis	loop	0/60
2 Bacillus subtilis	cylinder	0/60
2 Clostridium sporogenes	loop	0/60
2 Clostridium sporogenes	cylinder	0/60
3 Bacillus subtilis	loop	0/60
3 Bacillus subtilis	cylinder	0/60
3 Clostridium sporogenes	loop	0/60
3 Clostridium sporogenes	cylinder	0/60

Deactivation of Organisms by Minncare® Vapors

Estimated spore concentration of filter suspended 1 inch above	<u>B. subtilis</u> spores/filters
Control: water	
10 ³ spores	9.6 × 10 ³ CFU*
10 ⁶ spores	> 3 × 10 ³ CFU
Minncare® solution (1%)	
10 ³ spores	<1 CFU
10 ⁴ spores	<1 CFU
10 ⁵ spores	<1 CFU
10 ⁶ spores	<1 CFU

* Colony forming units

Minncare® Toxicity Assessment Summary

Toxicity summary	Specimen	Results	Reference
Oral toxicity	Male rats	LD ₅₀ = 2.43 (2.04-2.88) g/kg	Litchfield & Wilcoxon
	Femal rats	LD ₅₀ = 2.10 (1.92-2.30) g/kg	
Inhalation toxicity	Male rats	Established lethal concentration	
	Female rats	LC = 13,439 mg/cubic meter	
Skin sensitivity	Mice	No reaction	Burkhard's Test
	Humans	No visible effects	
Mucos membranes	New Zealand rabbits	Effects completely gone within 7 days	HH Draize
Dermatological sensitivity qualities	White Guinea pigs	No difference between control and test group	Klugman & Magnusson
Intravenous toxicity	Male rats	LD ₅₀ = 212 mg/kg	Litchfield & Wilcoxon

Minncare Dry Fog™ - A Best Way for Clean Room Disinfection.

⇒ An Integral Part of Modern Clean Room Disinfection Procedures

⇒ EPA Registered

⇒ Match FDA , PIC/S request.

⇒ Superior "Dry" Delivery System

⇒ Better Dispersion than "Wet" Systems

⇒ Effective for Very Large Areas: 20 ~ 1000 M³

⇒ Better for Critical Applications



- ⇒ Less Preparation Work
- ⇒ Less Corrosion
- ⇒ Less Residue

Droplet Size effect


⇒ The Minncare Dry Fog System Produces 7.5 μm Droplets

–Droplets Behave Like Gas Particles


- Bounce Off Surfaces

- Do Not Wet Surfaces


–Float in the Air Due to Brownian Motion.



The Phenomenon of Dry Fog



⇒ Small droplets bounce and do not burst upon collision.



⇒ Large droplets burst and make things wet.

Better Dispersion

⇒ Minncare Dry Fog Creates a Vapor State, It Disperses Throughout the Room More Completely Than a Liquid Mist.

⇒ Better Dispersion Means More Surface Contact.

⇒ Even Hard to Reach Surfaces, Like Behind Cabinets, Under Tables.

⇒ Less Condensation.

Relative Humidity

- ⇒ Initial Relative Humidity Should Be <50%.
 - ⇒ Optimum Activity Between 70 and 80% Relative Humidity.(1)
 - Measure Using [Optional Cable-Free Thermo Hygrometer](#).
 - ⇒ Maintain Relative Humidity Below 90% to Avoid Condensation
- (1) [Hideharu Shintani - Agents of Sterilization, Disinfection, and Antisepsis Used for Medical and Food Affiliation. National Institute of Health Sciences, Japan](#)

Room Size

- ⇒ The Number of Nozzles or Machines Required is Dependent on the Size of the Room
- ⇒ Each Nozzle Can Fog Up To 8,750 ft³ (250 m³) Effectively
- ⇒ Each Machine Can Fog Up To 35,000 ft³ (1,000 m³) With 4 Nozzles

Machine Positioning

- ⇒ Strategically Place the Machine in the Room
- ⇒ Aim Nozzles Away from Nearby Walls and Equipment
- ⇒ Position the Machines to Allow an Easy Flow of the Fog to All Parts of the Area

General Information

⇒ SETUP

- Amount of Minncare: 1.5 ml / M³ of Room Volume
- Total Solution Volume: Function of the Initial RH
- Ideal Initial Relative Humidity: less than 55%
- Ideal Room Temperature: 20-25°C

⇒ PROCESS

- Fog Dispersion Time: 15 - 120 min
- Exposure Time: 1 hour
- Venting Time: 0.5 - 2 hours
- ⇒ **TOTAL PROCESS TIME: 2 - 5 hours (including venting)**

※Notes: 1) All individuals MUST evacuate the room during the process.

Assure Safe Operation with the [Dry Fog Remote Control Unit](#)


2) After confirmation of chemical vapor residuals below limits ([using the Minncare Dry Fog Vapor Detection System](#)) the room may be returned to service.

The General Procedure

- Unpack & Assemble the System
- Calculate the Needed Water and Minncare Quantity
- Position Fogger
- Add Water
- Connect Pressurized Air Source
- Test System Operation

- Adjust Gauges
- Add Minncare
- Start Spraying

Dry Fog Setup – Software Support-Initial Calculations:



Position & Operate the System



Dry Fog – Air Residuals

⇒ Measure Using the [Minnicare Dry Fog Vapor Detection System](#)



⇒ Permissible Exposure Limits (PEL)

–Hydrogen Peroxide: 1 ppm

–Acetic Acid: 10 ppm

✳ Residuals on the surfaces : Much less than 0.00002 ml Acetic Acid / cm²



⇒ EXAMPLE - 70 M³



**6 log reduction in spore
B. Stearothermophilus
ATCC 7953**

**Surface control: 0 germ
Air control: 0 germ**

⇒ EXAMPLE - 950 m³



**6 log reduction in spore
B. Stearothermophilus
ATCC 7953**

**Surface control: 0 germ
Air control: 0 germ**

Evaluation

Customer Evaluation of Dry Fog

- ⇒ Diffusion Testing
- ⇒ Microbiological Testing

How Do Customers Validate Their Disinfection Process?

⇒ Tools

- Spores Strips: Paper, SS, Glass, Plastic
- Inoculated Petri Dishes
- Surface Microbiological Controls
- Air Microbiological Controls
- ✘ Cleanroom Disinfection Does Not Need a Specific Method or Specific Microorganisms for Process Validation. Controls Showing LRV Are Sufficient.

⇒ Possible Microorganisms for Cleanroom Disinfection Process Validation

- *Bacillus thuringiensis* (subtilis) ATCC 9372
- *Geobacillus stearothermophilus* ATCC 7953
- *Bacillus cereus* CIP783
- *Anterococcus hirae* CIP5855
- Use Caution with *B. stearothermophilus* ATCC 12980
- ✘ Some Versions Packaged for Autoclave Tests (Sealed from Vapor)

The Complete Solution for Clean Room Aerosol-Based Disinfection

MINNCARE DRY FOG SYSTEM

Current Methods of Room Disinfection

- ⇒ Surface Wiping
- ⇒ Manual Spraying (with Spray Bottles)
- ⇒ Heating Process (Vaporization)
- ⇒ Cold Process - Wet fogging
 - **Dry fogging**

✘ Disinfection is NOT a Cleaning Process

Currently Used Chemical Products


- ⇒ Alcohol based
- ⇒ Formaldehyde based
- ⇒ Gluteraldehyde based
- ⇒ Quaternary Ammonium Compounds
- ⇒ Cocktails of the above

Current Problems

- ⇒ Efficiency (more or less broad spectrum biocide)
- ⇒ Adaptation of the Micro-Organisms
- ⇒ Toxicity of the Chemicals Used
- ⇒ Long Contact Time
- ⇒ Very Long Time for Venting
- ⇒ Neutralization Needed
- ⇒ Corrosion
- ⇒ Variability of Chemical Application
- ⇒ Residuals



Minnicare® Cold Sterilant



⇒ An alternative solution for room

disinfection in Pharmaceutical,
Biotech, Medical Facilities,...

Benefits of Minncare®

- ⇒ Superior Biocidal Activity
- ⇒ All Components are Ultra-pure, Pharmaceutical Quality
- ⇒ US EPA (Environment Protection Agency) Registered Sterilant
- ⇒ No Heavy Metal Trace Contamination / Stabilizers.
- ⇒ Biodegradable – Decomposes to Acetic Acid, Water and Oxygen
- ⇒ No Toxic Aldehyde Type Vapors, Easy to Vent
- ⇒ Validated Residual Vapor Detection System

Minncare® Applications

- ⇒ Water Systems Disinfection
Tanks, Piping, Resins, Filters, RO Membranes,...
- ⇒ Surface Disinfection (Wiping and/or Manual Spraying)
- ⇒ Fogging



Minncare® Biocidal Activity

Minncare® is a potent antimicrobial agent and is effective against a broad spectrum of microorganisms, including:

- ⇒ Bacteria
- ⇒ Yeast and Molds
- ⇒ Mycobacteria
- ⇒ Bacteria Spores
- ⇒ Viruses

Activities of the most important biocides (Guyader, 1996)							
Biocides	Bacteria		Mycobacteria	Spores	Moulds	Yeasts	Virus
	Gram -	Gram +					
Peracetic acid	+++	+++		++	++	++	++
Alcohols	++	++		0	++	++	+
Alcohol (70?)	++	++	0	+	+	++	+
Glutaraldehyde	+++	+++	++	+	+++	++	++
Quat Ammonium	+++	+*	0	0	+	+	+
Chlorine	+++	+++	++	++	++	++	++
Hydrogen Peroxyde	+++	+++		+	+	+	0
Iodine	+++	+++	++	++	++	++	++

* Not active on *Pseudomonas*

Minncare® Mechanism of Action

- ⇒ Disrupts Sulfhydryl (-SH) and Sulfur (S-S) bonds in proteins and enzymes Important components in cells and membranes are broken by oxidative disruption.
- ⇒ Impede cellular activity by disrupting chemosmotic function of lipoprotein cytoplasmic membrane transport through rupture or dislocation of cell walls

- ⇒ Denature the properties of protein components by altering the nucleic acid structure of organisms.
- ⇒ Damage vegetative cells by oxidation with hydroxy radicals.
- ⇒ Produces organic radicals that act as reducing agents for spores.

“NO CHANCE” For Micro-Organisms to Build a Resistance

Minnicare® In Vitro Aqueous Testing

Species	Minnicare® Concentration .5%	Count per ml	Time for 100% kill in minutes
Bac. subtilis		6 × 10 ⁵	2.5
Bac: stearothermophilus		6 × 10 ⁶	2.5
Bac. subtilis NCTC 3610		2.4 × 10 ⁹	5.0
Bac. mesentericus		1.6 × 10 ⁹	5.0
Clostr. perfringens		1 × 10 ⁷	10.0
Clostr. tyrobutyricum		1 × 10 ⁷	5.0
Sacchar. cereisiae		6 × 10 ⁷	0.5
Cand. mycoderma		1.4 × 10 ⁸	0.5
Hansenula anomala		6.4 × 10 ⁸	0.5
Pichia membronaefaciens		4.8 × 10 ⁸	0.5
Pen. cameronense		1.7 × 10 ⁸	2.5
Mucor plumbeus		3 × 10 ⁶	2.5
Geotrichum candidum		2 × 10 ⁷	0.5
Byssoschlamys nivea		6 × 10 ⁷	0.5
Staph. aureus		2.6 × 10 ⁹	0.5
Strept. faecalis		4.6 × 10 ⁹	0.5
Kleb. aerogenes		2.3 × 10 ⁹	0.5
Ps. fluorescens		4.6 × 10 ⁹	0.5
Ps. aeruginosa		2 × 10 ⁹	0.5
Salm. thyphimurium		2.8 × 10 ⁹	0.5
Coryneb. rubrum		1 × 10 ⁷	1.0
Leuconostoc spec.		5.3 × 10 ⁹	0.5
Lactob. brevis		1.8 × 10 ⁹	0.5

Minnicare® AOAC Sporicidal Test

A total of 720 tests were performed on three lots of the test chemical. Results of the test are shown on the chart below.

Lot/Organism	Carrier type	Positives/Total
1 Bacillus subtilis	loop	0/60
1 Bacillus subtilis	cylinder	0/60
1 Clostridium sporogenes	loop	0/60
1 Clostridium sporogenes	cylinder	0/60
2 Bacillus subtilis	loop	0/60
2 Bacillus subtilis	cylinder	0/60
2 Clostridium sporogenes	loop	0/60
2 Clostridium sporogenes	cylinder	0/60
3 Bacillus subtilis	loop	0/60
3 Bacillus subtilis	cylinder	0/60
3 Clostridium sporogenes	loop	0/60
3 Clostridium sporogenes	cylinder	0/60

Deactivation of Organisms by Minncare® Vapors

Estimated spore concentration of filter suspended 1 inch above	<u>B. subtilis</u> spores/filters
Control: water	
10 ³ spores	9.6 × 10 ³ CFU*
10 ⁶ spores	> 3 × 10 ³ CFU
Minncare® solution (1%)	
10 ³ spores	<1 CFU
10 ⁴ spores	<1 CFU
10 ⁵ spores	<1 CFU
10 ⁶ spores	<1 CFU

* Colony forming units

Minncare® Toxicity Assessment Summary

Toxicity summary	Specimen	Results	Reference
Oral toxicity	Male rats	LD ₅₀ = 2.43 (2.04-2.88) g/kg	Litchfield & Wilcoxon
	Femal rats	LD ₅₀ = 2.10 (1.92-2.30) g/kg	
Inhalation toxicity	Male rats	Established lethal concentration	
	Female rats	LC = 13,439 mg/cubic meter	
Skin sensitivity	Mice	No reaction	Burkhard's Test
	Humans	No visible effects	
Mucos membranes	New Zealand rabbits	Effects completely gone within 7 days	HH Draize
Dermatological sensitivity qualities	White Guinea pigs	No difference between control and test group	Klugman & Magnusson
Intravenous toxicity	Male rats	LD ₅₀ = 212 mg/kg	Litchfield & Wilcoxon

Minncare Dry Fog™ - A Best Way for Clean Room Disinfection.

⇒ An Integral Part of Modern Clean Room Disinfection Procedures

⇒ EPA Registered

⇒ Match FDA , PIC/S request.

⇒ Superior "Dry" Delivery System

⇒ Better Dispersion than "Wet" Systems

⇒ Effective for Very Large Areas: 20 ~ 1000 M³

⇒ Better for Critical Applications



- ⇒ Less Preparation Work
- ⇒ Less Corrosion
- ⇒ Less Residue

Droplet Size effect


⇒ The Minncare Dry Fog System Produces 7.5 μm Droplets

–Droplets Behave Like Gas Particles


- Bounce Off Surfaces

- Do Not Wet Surfaces


–Float in the Air Due to Brownian Motion.



The Phenomenon of Dry Fog



⇒ Small droplets bounce and do not burst upon collision.



⇒ Large droplets burst and make things wet.

Better Dispersion

⇒ Minncare Dry Fog Creates a Vapor State, It Disperses Throughout the Room More Completely Than a Liquid Mist.

⇒ Better Dispersion Means More Surface Contact.

⇒ Even Hard to Reach Surfaces, Like Behind Cabinets, Under Tables.

⇒ Less Condensation.

Relative Humidity

- ⇒ Initial Relative Humidity Should Be <50%.
 - ⇒ Optimum Activity Between 70 and 80% Relative Humidity.(1)
 - Measure Using [Optional Cable-Free Thermo Hygrometer](#).
 - ⇒ Maintain Relative Humidity Below 90% to Avoid Condensation
- (1) [Hideharu Shintani - Agents of Sterilization, Disinfection, and Antisepsis Used for Medical and Food Affiliation. National Institute of Health Sciences, Japan](#)

Room Size

- ⇒ The Number of Nozzles or Machines Required is Dependent on the Size of the Room
- ⇒ Each Nozzle Can Fog Up To 8,750 ft³ (250 m³) Effectively
- ⇒ Each Machine Can Fog Up To 35,000 ft³ (1,000 m³) With 4 Nozzles

Machine Positioning

- ⇒ Strategically Place the Machine in the Room
- ⇒ Aim Nozzles Away from Nearby Walls and Equipment
- ⇒ Position the Machines to Allow an Easy Flow of the Fog to All Parts of the Area

General Information

⇒ SETUP

- Amount of Minncare: 1.5 ml / M³ of Room Volume
- Total Solution Volume: Function of the Initial RH
- Ideal Initial Relative Humidity: less than 55%
- Ideal Room Temperature: 20-25°C

⇒ PROCESS

- Fog Dispersion Time: 15 - 120 min
- Exposure Time: 1 hour
- Venting Time: 0.5 - 2 hours
- ⇒ **TOTAL PROCESS TIME: 2 - 5 hours (including venting)**

※Notes: 1) All individuals MUST evacuate the room during the process.

Assure Safe Operation with the [Dry Fog Remote Control Unit](#)


2) After confirmation of chemical vapor residuals below limits ([using the Minncare Dry Fog Vapor Detection System](#)) the room may be returned to service.

The General Procedure

- Unpack & Assemble the System
- Calculate the Needed Water and Minncare Quantity
- Position Fogger
- Add Water
- Connect Pressurized Air Source
- Test System Operation

- Adjust Gauges
- Add Minncare
- Start Spraying

Dry Fog Setup – Software Support-Initial Calculations:



Position & Operate the System



Dry Fog – Air Residuals

⇒ Measure Using the [MinnCare Dry Fog Vapor Detection System](#)



⇒ Permissible Exposure Limits (PEL)

–Hydrogen Peroxide: 1 ppm

–Acetic Acid: 10 ppm

✳ Residuals on the surfaces : Much less than 0.00002 ml Acetic Acid / cm²



⇒ EXAMPLE - 70 M³



**6 log reduction in spore
B. Stearothermophilus
ATCC 7953**

**Surface control: 0 germ
Air control: 0 germ**

⇒ EXAMPLE - 950 m³



**6 log reduction in spore
B. Stearothermophilus
ATCC 7953**

**Surface control: 0 germ
Air control: 0 germ**

Evaluation

Customer Evaluation of Dry Fog

- ⇒ Diffusion Testing
- ⇒ Microbiological Testing

How Do Customers Validate Their Disinfection Process?

⇒ Tools

- Spores Strips: Paper, SS, Glass, Plastic
- Inoculated Petri Dishes
- Surface Microbiological Controls
- Air Microbiological Controls
- ✘ Cleanroom Disinfection Does Not Need a Specific Method or Specific Microorganisms for Process Validation. Controls Showing LRV Are Sufficient.

⇒ Possible Microorganisms for Cleanroom Disinfection Process Validation

- *Bacillus thuringiensis* (subtilis) ATCC 9372
- *Geobacillus stearothermophilus* ATCC 7953
- *Bacillus cereus* CIP783
- *Anterococcus hirae* CIP5855
- Use Caution with *B. stearothermophilus* ATCC 12980
- ✘ Some Versions Packaged for Autoclave Tests (Sealed from Vapor)